

# **AZtecCrystal**

## A new era in EBSD data processing

Electron backscatter diffraction (EBSD) is an increasingly essential microanalysis technique, driven by high-speed CMOS-based detectors and a growth in markets in which EBSD helps to solve real materials problems. The result is a requirement to increase the breadth and accessibility of data interrogation tools.

AZtecCrystal sets a new standard for EBSD data analysis: a totally new processing package, designed for the EBSD novice and expert alike, that will unlock the hidden characteristics in every sample.

### Key features include:

- Single intuitive yet flexible user interface
- Seamless integration with the AZtec acquisition platform
- Newly developed HDF5 data format including both EDS and EBSD results
- Full map analysis with multiple customisable displays
- Full texture analysis using pole figures and orientation distribution functions
- Multiple advanced tools, including materials properties, machine-learning based classification, data subsetting, comprehensive grain sizing and disorientation analysis

#### **Flexible**

- Multi modal display: view the analytical tools relevant to every materials problem
- Multi layered map creation with no component limits
- Ideal for all computers, from laptops toworkstations

#### **Fast**


- Quaternion based algorithms for fastest data processing
- Handles any size dataset from AZtec with ease: even the largest datasets are processed within seconds, boosted by the batch processing tool
- Relevant data and settings are stored within each project for fast reloading

#### Intelligent

- Designed around intelligent routines, such as void detection, GND density calculation, parent grain size determination and disorientation angle congruence
- Comprehensive grain size analysis to ASTM standards with over 10 display options
- Smart defaults and settings for all materials

#### **New Science**

- AZtecCrystal makes advanced EBSD accessible to all users, opening up new research opportunities
- Dynamic platform multiple new tools and features to be incorporated in future releases



The materials presented here are summary in nature, subject to change, and intended for general information only. Additional details are available. Oxford Instruments NanoAnalysis is certified to ISO9001. ISO14001 and OHSAS 18001. AZtec and Symmetry are Registered Trademarks of Oxford Instruments plc, all other trademarks acknowledged. © Oxford Instruments plc, 2020. All rights reserved.

